Organic Initiative Ltd
PO Box 50165 Albany
Auckland 0752, New Zealand
Email: info@oi4me.com
US toll-free Ph: +64 9 212 6744
Visit oi4me.com for more information.

ENGLISH

The menstrual cup is an alternative to tampons and pads. Its function is to collect the flow of your period and its application is very similar to a tampon.

The Oi Cup™ is inserted a bit like a tampon and then pops open to form a cup which catches your menstrual flow. It can stay in for 6-12 hours. You then wash it and reuse it for up to three years.

The Oi Cup™ is available in three sizes to suit different women and their flow. The following is a guideline, but it may take time to work out what size fits best for your unique shape.

- **Small:** Suitable for girls and women under 18 years old, who do not have sex regularly, or women who are smaller inside.
- **Medium:** The most common size and suitable for women over 18 years old, and women who have not given birth vaginally.
- **Large:** For women who are larger inside, this may be because they have given birth vaginally.

USING THE OI CUP™

Getting used to a menstrual cup can be a little difficult at first, but be patient. After a little practice inserting and removing it, it is very easy and very comfortable. Oi Cup™ users love how simply and effectively it works.

HOW TO USE THE OI CUP™

1. Sterilize the cup

Before using the cup for the first time each month you must sterilize it:

- First wash your hands.
- Bring water to the boil in a saucepan and dip the cup in the water for 3-4 minutes, covering it completely.
- Remove the cup being careful not to burn yourself, and it will be ready to use once it has cooled to your skin temperature.

2. Get yourself relaxed and ready

- Wash your hands thoroughly.
- Adopt a comfortable position. You can sit on the toilet, squat or stand.
- Fold the cup to make it easier to insert:

The simplest option to the fold the cup is:

C Fold

Press both sides of the cup to put them together and then fold the cup in half forming a 'C'.


PRESS
ES: PRENSA
FR: PRESSE

FOLD
ES: DOBLEZ
FR: PLIER

3. Insert the Oi Cup™


- While holding the cup folded with one hand, with the other hand separate your labial lips and push the cup into your vagina and towards the tailbone or base of your spine. The cup sits a bit lower than a tampon, but not so low that the stem becomes uncomfortable.

INSERT
ES: INTRODUCIR
FR: INSÉREZ

- Once you have inserted the cup, let it loose, it will 'pop' open and perfectly fit your vaginal wall.

- To be sure that the cup has created a seal, turn it slightly by grabbing the base of the cup just above the stem.

TWIST
ES: GIRO
FR: TORSION

4. Removing it to empty

- You must empty the cup at least once every 12 hours for your personal health.
- Once you are familiar with its use and your flow, you will know how frequently you need to empty it.
- If your flow is heavier, you may need to empty it sooner. Typically a cup is emptied after six hours and before 12 hours after placement.

How to remove it

- Get into a comfortable position (sitting on the toilet, squatting or standing).
- Place your index finger and thumb in the vagina to locate the stem, reach just above the stem, pinch the bottom of the cup to release suction and gently pull down and out keeping it in an upright position.
- Push down with your abdominal/pelvic muscles to make it easier.
- Gently pull the stem until you are able to firmly hold the base of the cup.
- Remove the cup completely, while still holding it vertically so that the content remains within it.
- Empty the content into the toilet, or the shower if this is where you prefer to remove it.

Re-insert

- Wash the cup with clean water, preferably warm, so that it is ready for using again. Remember to thoroughly wash the suction holes.
- You can also use neutral pH soap without perfume. If you do not

have access to clean water you can also use toilet paper or a handkerchief to dry and then wash it with clean water later.

Sterilize and store

- After your period is finished, we recommend sterilizing again to keep it perfectly clean.
- Store in the Oi cotton canvas bag provided.
- Next time you use it, you need to sterilize it again before inserting.

TIPS

1. Avoid discomfort: Your cup should be very comfortable. If it isn't then it may not be positioned correctly.

- Remove the cup and try the process again with a slightly different angle.
- You may find the stem too long, generally it shouldn't protrude outside the vagina. It is OK to trim it until you find the most comfortable length for you, as long as you don't damage the base of the cup.

2. Removing the cup: This can be tricky at first, but don't panic if it seems a bit hard, try these tips:

- Relax and stay calm.
- Breathe in and out slowly.
- Push your vagina/pelvic muscles several times so that if the cup is too high in the vagina it falls through.
- Try holding the base of the cup rather than the stem and making small movements from side to side as you roll it gently out.
- Press the cup base with your thumb and index finger for air to enter and help release the suction.
- **No leaks: Your cup should not leak. If yours is leaking it is probably not placed in correctly, or it is full.**
 - It should be placed at the entrance of the vagina, lower than a tampon.
- You should be able to rotate it easily inside you.
- Check the cup size, perhaps you have not chosen the right size and you need a larger one from our Oi Cup™ range.

- If none of these tips help, and you use a large size, you may need some exercises to strengthen the pelvic floor muscles. Consult with your Doctor. Or switch to a tampon or pad.

ANSWERS TO FREQUENTLY ASKED QUESTIONS

- **NO.** It won't disappear inside you.
- **NO.** You don't need to remove the cup to urinate or make bowel movements.
- **NO.** Your cup shouldn't leak if you have the right size and have placed it correctly.
- **YES.** You can use the cup at night for up to 12 hours.
- **YES.** The cup can be used while you exercise and play sport including swimming.
- **YES.** The cup is compatible with the IUD or diaphragm.
- **NEVER.** Use the cup as a contraceptive device, it isn't one.
- **SEEK.** Medical attention if you have any questions about TSS or cup use.

WARNINGS

- Use the cup according to the instructions in this leaflet and only during your period.
- The Oi Cup™ is a product of intimate use so you should not share it.
- For cleaning do not use abrasive cleaners or disinfectants such as tea tree oil, Aleppo or Alep soap.
- Do not use oil-based lubricants to insert the cup into the vagina.
- Never keep the cup in a plastic bag or airtight containers.
- If you travel, use safe drinking water for washing.
- In case of burning, irritation, genital or vaginal inflammation, discomfort during urination or other symptoms that you think are abnormal, immediately remove the cup and ask your doctor.
- Symptoms that warn of possible TSS are: sudden high fever (39 °C/102 °F or more), dizziness, near fainting when standing up, diarrhea, vomiting, fainting, muscle aches, and skin irritation similar to sunburn. It may be that not all symptoms are present, only a few. In more advanced stages of the disease, the skin may begin to peel. TSS can occur during menstruation or shortly thereafter. It can move quickly with symptoms similar to flu, even to a serious fatal disease.
- If you have any of these symptoms while using the menstrual cup, remove it and talk with your doctor immediately. If you already have had TSS do not use the cup.
- Sterilize the cup again before inserting.

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